



### Resin Cements for Use with Empress and Eris

Brand	Manufacturer	Empress 1 Staining	Empress 1 Layering	Empress Esthetic Veneer	Eris for E2
RelyX ARC	3M ESPE	Yes	Yes	Yes	Yes
RelyX Veneer	3M ESPE	Yes	Yes	Yes	Yes
C&B Cement	Bisco	Yes	Yes	Yes	Yes
Choice	Bisco	Yes	Yes	Yes	Yes
Illusion	Bisco	Yes	Yes	Yes	Yes
Insure	Cosmedent	Yes	Yes	Yes	Yes
Prevue	Cosmedent	Yes	Yes	Yes	Yes
Ultra-Bond Quik	Denmat	Yes	Yes	Yes	Yes
Calibra	Dentsply	Yes	Yes	Yes	Yes
Appeal	Ivoclar Vivadent	Yes	Yes	Yes	Yes
Multilink	Ivoclar Vivadent	Yes	Yes	Yes	Yes
Variolink II	Ivoclar Vivadent	Yes	Yes	Yes	Yes
Nexus	Kerr	Yes	Yes	Yes	Yes
Panavia 21	Kuraray	Yes	Yes	Yes	Yes
Panavia F	Kuraray	Yes	Yes	Yes	Yes
C&B Metabond	Parkell	Yes	Yes	Yes	Yes
Cement-it	Pentron	Yes	Yes	Yes	Yes
Lute-it	Pentron	Yes	Yes	Yes	Yes

### Conventional Cements for Use with Eris for E2

VivaGlass CEM	Ivoclar Vivadent	No	No	No	Yes
Fuji 1	GC	No	No	No	Yes
Ketac CEM	3M ESPE	No	No	No	Yes
Relyx Unicem	3M ESPE	No	No	No	Yes

**Resin Reinforced or Resin Modified Hybrid Glass Ionomer Cements should not be used with IPS Empress or Eris for E2 due to expansion**

RelyX Luting Cement	3M ESPE	No	No	No	No
Fuji CEM	GC	No	No	No	No
Fuji Plus	GC	No	No	No	No

**Empress Cements Prepared by Ivoclar Vivadent 2/21/2005**